



CLINICAL AND
LABORATORY
STANDARDS
INSTITUTE

2nd Edition

CLSI C38™

Control of Preexamination Variation in Trace Element Determinations

CLSI C38 provides guidelines for patient preparation, specimen collection, transport, and processing for the measurement of trace elements in a variety of biological matrixes.

A guideline for global application developed through the Clinical and Laboratory Standards Institute consensus process.

Clinical and Laboratory Standards Institute

Setting the standard for quality in medical laboratory testing around the world.

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advances in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

Appeal Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeal, documented in the CLSI *Standards Development Policies and Processes*, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute

P: +1.610.688.0100

F: +1.610.688.0700

www.clsi.org

standard@clsi.org

Control of Preexamination Variation in Trace Element Determinations

Robert L. Jones, PhD
Lisa Strong, PhD
Sandra Fitzpatrick, MAACB
Mitch André Garcia, PhD
Jeffery M. Jarrett, MS
Kamisha Johnson-Davis, PhD, DABCC, FADLM
Magen Martinez

Fred Leland McClure, MSci, PhD, F-ABFT
Michelle McLean, MS, MLS(ASCP), BS
Riley Murphy, PhD
Patrick J. Parsons, PhD, CChem, FRSC
Denise Rokke
Carl E. Wolf, PhD, MS, F-ABFT

Abstract

Clinical and Laboratory Standards Institute C38—*Control of Preexamination Variation in Trace Element Determinations* is intended for persons responsible for the collection and processing of samples used for trace, toxic, and essential metal and metalloid determinations. The guideline covers patient preparation, as well as considerations for collection, transport, and processing of specimens by element. Contamination control and quality assurance programs are also discussed.

Clinical and Laboratory Standards Institute (CLSI). *Control of Preexamination Variation in Trace Element Determinations*. 2nd ed. CLSI guideline C38 (ISBN 978-1-68440-240-3 [Print]; ISBN 978-1-68440-241-0 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2024.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at www.clsi.org.

If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:

P: +1.610.688.0100 **F:** +1.610.688.0700 **E:** customerservice@clsi.org **W:** www.clsi.org

Copyright ©2024 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, or other product or material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to permissions@clsi.org.

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedures manual at a single site. To request permission to use this publication in any other manner, e-mail permissions@clsi.org.

To read CLSI's full Copyright Policy, please visit our website at <https://clsi.org/terms-of-use/>.

Suggested Citation

CLSI. *Control of Preexamination Variation in Trace Element Determinations*. 2nd ed. CLSI guideline C38. Clinical and Laboratory Standards Institute; 2024.

Previous Edition:
September 1997

CLSI C38-Ed2

ISBN 978-1-68440-240-3 (Print)

ISBN 978-1-68440-241-0 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 44, Number 22

.....

Committee Membership

Consensus Council

The Consensus Council sets priorities for CLSI standards development and votes on Final Draft documents to confirm that process requirements have been met. Consensus Council members are listed on the CLSI website: <https://clsi.org/standards-development/consensus-council/>

Document Development Committee on Control of Preanalytical Variation in Trace Element Determinations

Robert L. Jones, PhD
Chairholder
USA

Kamisha Johnson-Davis, PhD, DABCC,
 FADLM
 ARUP Laboratories/University of Utah
 USA

Denise Rokke
 Mayo Clinic
 USA

Lisa Strong, PhD
Vice-Chairholder
Minnesota Department of Health
USA

Magen Martinez
 Tricore Reference Laboratories
 USA

Carl E. Wolf, PhD, MS, F-ABFT
 Virginia Commonwealth University
 Health System
 USA

Sandra Fitzpatrick, MAACB
 LabPlus Auckland District Health Board
 New Zealand

Michelle McLean, MS, MLS(ASCP), BS
 Roche Diagnostics
 USA

Mitch André Garcia, PhD
 BD
 USA

Patrick J. Parsons, PhD, CChem, FRSC
 New York State Department of Health
 USA

Expert Panel on Chemistry and Toxicology

Expert panel volunteers support the development of CLSI documents by providing technical expertise in specialty areas. Expert panel members are listed by area of expertise on the CLSI website: <https://clsi.org/standards-development/expert-panels/>

Acknowledgment

CLSI, the Consensus Council, and the Document Development Committee on Control of Preanalytical Variation in Trace Element Determinations gratefully acknowledge the following volunteers for their important contributions to the revision of this guideline:

Jeffery M. Jarrett, MS
Centers for Disease Control and
Prevention
USA

Fred Leland McClure, MSci, PhD,
F-ABFT
Drug Testing Consultations
USA

Riley Murphy, PhD
NMS Labs
USA

Contents

Abstract	i
Committee Membership	iii
Foreword	vii
Chapter 1: Introduction	1
1.1 Scope	2
1.2 Background	2
1.3 Standard Precautions	2
1.4 Terminology	3
Chapter 2: Path of Workflow	7
2.1 Contamination Control	8
2.2 Specimen Selection	12
Chapter 3: Specific Elements	19
3.1 Aluminum	20
3.2 Antimony	22
3.3 Arsenic	24
3.4 Barium	25
3.5 Bismuth	26
3.6 Cadmium	27
3.7 Chromium	29
3.8 Cobalt	31
3.9 Copper	33
3.10 Gadolinium	34
3.11 Gold	35
3.12 Iodine	36
3.13 Iron	38
3.14 Lead	39
3.15 Manganese	42
3.16 Mercury	43
3.17 Molybdenum	45
3.18 Nickel	46
3.19 Selenium	47
3.20 Titanium	48
3.21 Uranium	49

Contents (Continued)

3.22 Vanadium.....	51
3.23 Zinc.....	52
Chapter 4: Elemental Speciation and Fractionation	55
4.1 Preexamination Conditions for Arsenic Speciation Analytical Methods	56
4.2 Preexamination Conditions for Mercury Speciation Analytical Methods	58
Chapter 5: Analytical Instrumentation.....	61
5.1 Description of Instrumentation.....	62
5.2 Calibration	63
5.3 Verifying the Required Measurement Accuracy at Defined Concentrations	64
Chapter 6: Conclusion.....	65
Chapter 7: Supplemental Information	67
References	68
Additional Resources	78
Appendix A. Summary of Elements.....	79
Appendix B. Proficiency Testing Programs	80
Appendix C. Recommended Reporting Units.....	82
Appendix D. Biological Reference Interval Studies.....	83
The Quality Management System Approach.....	98

Foreword

Preexamination factors are probably the most important cause of erroneous trace element reference data in biological matrixes today. Developing sensitive, specific, and accurate analytical technology at an acceptable cost has moved the determination of trace, toxic, and essential metals from research facilities into a wide range of medical laboratories. Expanding knowledge of trace element nutrition and toxicity has increased clinical demand for these assays. However, with increased sensitivity and lower limits of detection, the problem of specimen contamination with the metal of interest has been magnified greatly. It is vital that the accurately determined trace, toxic, or essential element concentration reflects the condition of the patient and not contamination introduced during collection, processing, handling, and analysis.

Earlier attempts to define reference interval data for many of the trace, toxic, and essential elements provided intervals that were far wider than what are now accepted as a population “normal.” This resulted from a lack of awareness that the ubiquity of many trace elements in the environment required special precautions from preexamination processes through the actual analysis.

In this guideline, the components of specimen collection and preexamination and examination processing that can contribute to trace element contamination are covered. The trace elements most commonly tested for clinical purposes are individually listed. For each element, the optimal specimen for assessment, preexamination factors to consider in patient preparation and population reference intervals, or concentrations suggesting elevated exposures or deficiency, are described.

Overview of Changes

This guideline replaces CLSI C38-A, published in 1997. Several changes were made in this edition, including:

- Adding information on indications, patient preparations, sample selection, transport and storage, reference intervals, and conversion factors for each element
- Updating population reference intervals from the latest US National Health and Nutrition Examination Survey for most elements
- Adding a chapter on analytical instrumentation, recommended reporting units (see Appendix C) and an extensive appendix on reference interval studies (see Appendix D)

NOTE: The content of this guideline is supported by the CLSI consensus process and does not necessarily reflect the views of any single individual or organization.

KEY WORDS

analytical instrumentation

contamination control

essential element

metalloid

specimen collection

trace element

toxic element

Chapter ①

Introduction

Control of Preexamination Variation in Trace Element Determinations

1 Introduction

1.1 Scope

C38 provides guidance for preexamination considerations for trace element analysis in biological matrixes (ie, body fluids, such as whole blood, urine, hair, nails, human milk, and tissues). This guidance includes applicable recommendations for patient preparation (such as diurnal variation, diet, medications including chelation therapy and mineral supplements; the use of alcohol and tobacco; time of collection relative to the above activities, including workplace schedules, and homeostasis), specimen collection (most appropriate specimen matrix for screening, diagnosis, or monitoring), use of anticoagulants and/or preservatives; good practices for collection procedures (including containers and devices), transport requirements, contamination monitoring and control (including water, reagents, and examination consumables), and specimen processing. Specific reference is made to those elements that are known to be essential or toxic for humans and are therefore likely to be measured for clinical reasons.

This guideline will benefit those in the medical, public health, and environmental laboratories as well as manufacturers and regulatory and accreditation agencies.

This guideline will not include an extensive review of examination protocols and postexamination considerations, lead (Pb) testing (see CLSI C40¹) or content that is not consistent with current practices and guidelines.

1.2 Background

It is recognized that much of the pioneering research published in trace element literature is based on erroneously derived reference interval data.² The source of the problem was the lack of recognition of exogenous specimen contamination, which could have occurred at the collection, transport, processing, or examination stages. Thus, reference intervals for ultratrace elements, such as chromium (Cr), or acceptable blood concentrations for toxic elements, such as aluminum (Al), have decreased severalfold over the past two decades.

The use of increasingly sensitive methods, such as electrothermal atomic absorption spectrometry (ETAAS) or inductively coupled plasma mass spectrometry (ICP-MS); increasing interest in ultratrace elements; and the need for precise and accurate analyses for elements such as Pb, at extremely low levels, have accentuated the problems of examination and preexamination contamination.³

The intent of this guideline is to develop an awareness of the factors that affect the determination of trace elements in a variety of specimen types, foster communication between the laboratorian performing the test and those responsible for collecting the specimen, and provide suggested protocols for eliminating preexamination variability.

1.3 Standard Precautions

Because it is often impossible to know what isolates or specimens might be infectious, all patient and laboratory specimens are treated as infectious and handled according to “standard precautions.” Standard precautions are guidelines that combine the major features of “universal precautions and body substance isolation” practices. Standard precautions cover the transmission of all known infectious agents and thus are more comprehensive than universal precautions, which are intended to apply only to transmission of bloodborne pathogens. Published