

# EP25

## Evaluation of Stability of *In Vitro* Medical Laboratory Test Reagents

This guideline provides recommendations for establishing and verifying shelf-life and in-use stability claims for *in vitro* diagnostic medical laboratory test reagents such as reagent kits, calibrators, and control products.

.....

A guideline for global application developed through the Clinical and Laboratory Standards Institute consensus process.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advances in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeal Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeal, documented in the CLSI *Standards Development Policies and Processes*, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute

P: +1.610.688.0100

F: +1.610.688.0700

[www.clsi.org](http://www.clsi.org)

[standard@clsi.org](mailto:standard@clsi.org)

---

# Evaluation of Stability of *In Vitro* Medical Laboratory Test Reagents

Karl De Vore, BA, SSBB  
Mark Holland, PhD  
Lorin M. Bachmann, PhD, DABCC, MT(ASCP)  
John G. Chromczak, MLS(ASCP)SH  
Florence Crepet  
Julian Duncan, MA, PhD  
Sheryl Fennell, MT(ASCP)  
Justine Fertey

Qing H. Meng, MD, PhD, DABCC, FAACC  
James F. Pierson-Perry  
Michael Pistorino, MS  
Brian Renley, MS  
Vickie Lynn Simmons, BS  
Nadeem A. Tusneem, MA  
Anca Roxana Varlan, PhD  
Kong M. Xiong, PhD

## Abstract

---

Clinical and Laboratory Standards Institute guideline EP25—*Evaluation of Stability of In Vitro Medical Laboratory Test Reagents* provides recommendations and regression-based procedures for establishing and subsequently confirming stability-related claims of *in vitro* medical laboratory reagents such as reagent kits, calibrators, control products, and sample diluents. This guideline was written primarily for manufacturers and regulatory agencies but will also be of interest to medical laboratories and developers of laboratory-developed test methods. It provides information on the design, implementation, data analysis, and documentation needs for studies to establish and confirm shelf life and in-use life of *in vitro* diagnostic products. Additional topics cover the assessment of product transport conditions on stability, use of mean kinetic temperature to reflect product exposure to temperature changes during distribution and storage, and accelerated stability testing.

Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Stability of In Vitro Medical Laboratory Test Reagents*. 2nd ed. CLSI guideline EP25 (ISBN 978-1-68440-178-9 [Print]; ISBN 978-1-68440-179-6 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2023.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org).

**If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:**

**P:** +1.610.688.0100 **F:** +1.610.688.0700 **E:** [customerservice@cls.org](mailto:customerservice@cls.org) **W:** [www.clsi.org](http://www.clsi.org)

Copyright ©2023 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, derivative product, or other material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedures manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

## Suggested Citation

CLSI. *Evaluation of Stability of In Vitro Medical Laboratory Test Reagents*. 2nd ed. CLSI guideline EP25. Clinical and Laboratory Standards Institute; 2023.

**Previous Edition:**

September 2009

EP25-Ed2

ISBN 978-1-68440-178-9 (Print)

ISBN 978-1-68440-179-6 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 43, Number 8

.....

## Committee Membership

### Consensus Council

The Consensus Council sets priorities for CLSI standards development and votes on Final Draft documents to confirm that process requirements have been met. Consensus Council members are listed on the CLSI website: <https://clsi.org/standards-development/consensus-council/>

### Document Development Committee on Evaluation of Reagent Stability

**Karl De Vore, BA, SSBB**

**Chairholder**

**Bio-Rad Laboratories, Inc.**

**USA**

Julian Duncan, MA, PhD

Representing the World Health

Organization

England

Ravinder Jit Singh, PhD

Mayo Clinic

USA

**Mark Holland, PhD**

**Vice-Chairholder**

**Beckman Coulter**

**USA**

Qing H. Meng, MD, PhD, DABCC,

FAACC

The University of Texas MD Anderson

Cancer Center

USA

Anca Roxana Varlan, PhD

Questat, Inc.

Canada

Lorin M. Bachmann, PhD, DABCC,

MT(ASCP)

Virginia Commonwealth University

Health System

USA

Marion Niessner, PhD

Roche Diagnostics GmbH

Germany

Kathleen B. Whitaker, PhD

FDA Center for Devices and

Radiological Health

USA

Florence Crepet

Diagnostica Stago

France

James F. Pierson-Perry

USA

### Expert Panel on Evaluation Protocols

Expert panel volunteers support the development of CLSI documents by providing technical expertise in specialty areas. Expert panel members are listed by area of expertise on the CLSI website: <https://clsi.org/standards-development/expert-panels/>

### Staff

Clinical and Laboratory Standards

Institute

USA

Laura Martin

*Editorial Manager*

Kristy L. Leirer, MS

*Editor*

Tabitha Kern, MS, MLS(ASCP)<sup>CM</sup>

*Program Manager*

Catherine E.M. Jenkins, ELS

*Editor*

Lisa M.W. Walker, MS, ELS

*Editor*

## Acknowledgment

---

CLSI, the Consensus Council, and the Document Development Committee on Evaluation of Reagent Stability gratefully acknowledge the following volunteers for their important contributions to the revision of this guideline:

John G. Chromczak, MLS(ASCP)SH  
Diagnostica Stago  
USA

Michael Pistorino, MS  
Instrumentation Laboratory  
USA

Nadeem A. Tusneem, MA  
Q<sup>2</sup> Solutions – EA Genomics  
USA

Sheryl Fennell, MT(ASCP)  
Siemens Healthcare Diagnostics, Inc.  
USA

Brian Renley, MS  
Abbott  
USA

Kong M. Xiong, PhD  
Exact Sciences  
USA

# Contents

Abstract .....	i
Committee Membership .....	iii
Foreword .....	vii
<b>Chapter 1: Introduction .....</b>	<b>1</b>
1.1 Scope .....	2
1.2 Standard Precautions .....	2
1.3 Terminology .....	2
<b>Chapter 2: Overview of the Stability Testing Process .....</b>	<b>9</b>
2.1 Process Flow Chart .....	10
2.2 Overview of Stability Testing .....	11
2.3 Definition of Stability for <i>In Vitro</i> Diagnostic Products .....	12
2.4 Types of Stability Studies .....	14
2.5 Extending Product Shelf-Life Claims .....	15
2.6 Product Design Changes .....	15
2.7 Considerations for Products With Multiple Uses .....	16
2.8 Stability Testing for Qualitative Examinations .....	16
<b>Chapter 3: Stability Validation: Developing a Stability Study Design .....</b>	<b>17</b>
3.1 Stability Plan Overview .....	18
3.2 Elements of a Stability Plan .....	18
3.3 Transport Simulation Studies Plan .....	25
3.4 Use of Mean Kinetic Temperature .....	26
3.5 Combining Studies .....	28
<b>Chapter 4: Stability Validation: Implementing a Stability Study .....</b>	<b>29</b>
4.1 Stability Testing Experimental Steps .....	30
4.2 Data Quality Monitoring Investigations .....	30
<b>Chapter 5: Stability Validation: Analyzing Stability Study Data .....</b>	<b>33</b>
5.1 Technical Data Review .....	34
5.2 Data Analysis Steps for Each Metric's Dataset .....	35
5.3 Overall Stability Assessment .....	37
5.4 Study Conclusions and Claims .....	37
5.5 Stability Report .....	38

## Contents (Continued)

<b>Chapter 6: Verifying (Confirming) Established Stability Claims</b> .....	<b>39</b>
6.1 Rationales for Verifying Stability Claims .....	40
6.2 Preexamination Points to Consider .....	40
6.3 Examination Points to Consider .....	41
6.4 Postexamination Points to Consider .....	42
<b>Chapter 7: Special Topics</b> .....	<b>45</b>
7.1 Accelerated Stability Testing .....	46
7.2 Difficult Samples .....	61
7.3 Nontest Material Changeover Studies .....	62
<b>Chapter 8: Conclusion</b> .....	<b>63</b>
<b>Chapter 9: Supplemental Information</b> .....	<b>65</b>
<b>References</b> .....	66
<b>Appendix A.</b> Power Analysis for Stability Studies Based on Linear Regression .....	68
<b>Appendix B.</b> Example of Use of Arrhenius Equation With Accelerated Stability Testing Data to Predict Shelf Life of an <i>In Vitro</i> Diagnostic Control Product .....	76
<b>Appendix C.</b> Stability Study Examples .....	79
<b>Appendix D.</b> Statistical Basis for Stability Data Analysis .....	106
<b>The Quality Management System Approach</b> .....	112

## Foreword

Stability of an *in vitro* diagnostic (IVD) product reflects its ability to maintain consistent performance characteristics over time. Unlike precision, bias, and other common performance attributes, product stability is rarely assessed directly by end-user testing. As such, there is increased burden on manufacturers and developers of laboratory-developed tests (LDTs) to ensure that stability claims are developed from experimental designs and data analyses that are appropriate for each product's needs and applications.

Products, in the context of this guideline, represent end-user consumable products sold for performing laboratory measurements on patient samples or other samples claimed as appropriate by the manufacturer, such as veterinary or contrived matrixes made from biological materials. Examples of such products are IVD reagent kits or LDTs and their associated calibrators, controls, sample diluents, and system generic reagents. The information also applies to qualitative and semiquantitative or semiquantitative tests.

The content of this guideline is aligned with international standards for stability and internationally recognized guidance documents relative to stability study design and analyses.<sup>1-3</sup> Although these guidance documents were developed for drugs and drug substances, much of their content is directly relevant to IVD reagents.

### Overview of Changes

This guideline replaces the previous edition of the approved guideline, EP25-A, published in 2009. Several changes were made in this edition, including:

- Revising the approach to statistical power analysis for planning studies to assume there will be some drift in reagent performance
- Designing a testing plan to demonstrate the drift is within the allowable drift limit
- Eliminating the custom of using the t-test of regression slope results ( $P > 0.05$ ) as a rationale for passing a stability assessment. This practice tends to reward the manufacturer for designing underpowered stability studies.
- Eliminating the requirement for a confidence interval within the acceptance criteria as a basis for stating claims
- Revising the regression analysis approach that now requires point data at the claimed time and beyond the claimed time
- Expanding the practices for transport simulation stability testing
- Adding a discussion about the use of mean kinetic temperature as an integrated measure of temperature changes experienced by a product during distribution and storage
- Expanding the use and practices for accelerated stability testing

**NOTE:** The content of this guideline is supported by the CLSI consensus process and does not necessarily reflect the views of any single individual or organization.

#### KEY WORDS

**accelerated stability**

**in-use life**

**stability monitoring**

**allowable drift**

**isochronous design**

**stability plan**

**equivalence testing**

**mean kinetic temperature**

**transport simulation**

**expiry dating**

**shelf life**

This page is intentionally left blank.

# Chapter ①

## Introduction

# Evaluation of Stability of *In Vitro* Medical Laboratory Test Reagents

## 1 Introduction

### 1.1 Scope

This guideline provides information for establishing and verifying (confirming) shelf-life and in-use stability claims for quantitative *in vitro* diagnostic (IVD) medical laboratory reagents or products. The information also applies to qualitative IVD products, provided that an underlying continuous response or signal responsible for the qualitative result(s) is available to the investigator. This guideline also includes background information and typical content for creating a stability testing plan, determining the logistics for performing the studies, conducting recommended data analyses, and documenting stability claims for a product. Additional topics include assessment of product distribution conditions on stability claims (transport simulation), verification of stability claims, appropriate uses of accelerated testing, and considerations for testing with difficult samples.

The intended users of this guideline are primarily product manufacturers and regulatory agencies. Medical laboratorians may find this information useful for interpreting and confirming commercial product stability claims (eg, in-use life of QC materials), as well as for establishing stability attributes for laboratory-developed test methods. For this guideline, “products” is understood to encompass reagents, calibrators, controls, diluents, and similar materials that are used as IVD medical devices to conduct a measurement procedure for a measurand of medical interest.

This guideline does not cover instrument systems, laboratory equipment, software, or patient specimens. Stability testing of raw materials or components of reagent kits or consumables is not covered explicitly; however, the principles described in this guideline can be adapted by manufacturers for that purpose.

### 1.2 Standard Precautions

Because it is often impossible to know what isolates or specimens might be infectious, all patient and laboratory specimens are treated as infectious and handled according to “standard precautions.” Standard precautions are guidelines that combine the major features of “universal precautions and body substance isolation” practices. Standard precautions cover the transmission of all known infectious agents and thus are more comprehensive than universal precautions, which are intended to apply only to transmission of bloodborne pathogens. Published guidelines are available that discuss the daily operations of diagnostic medicine in humans and animals while encouraging a culture of safety in the laboratory.<sup>4</sup> For specific precautions for preventing the laboratory transmission of all known infectious agents from laboratory instruments and materials and for recommendations for the management of exposure to all known infectious diseases, refer to CLSI document M29.<sup>5</sup>

### 1.3 Terminology

CLSI, as a global leader in standardization, is firmly committed to achieving global harmonization whenever possible. Harmonization is a process of recognizing, understanding, and explaining differences while taking steps to achieve worldwide uniformity. CLSI recognizes that medical conventions in the global metrological community have evolved differently in different countries and regions and that legally required use of terms, regional usage, and different consensus timelines are all important considerations in the harmonization process. CLSI recognizes its important role in these efforts, and its consensus process focuses on harmonization of terms to facilitate the global application of standards and guidelines. Table 1 is provided to clarify the intended interpretations of the following terms.