

# H47

## One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test

This document provides guidelines for performing the PT and APTT tests in the clinical laboratory, for reporting results, and for identifying sources of error.

.....  
A guideline for global application developed through the Clinical and Laboratory Standards Institute consensus process.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advances in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeal Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeal, documented in the CLSI *Standards Development Policies and Processes*, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute

P: +1.610.688.0100

F: +1.610.688.0700

[www.clsi.org](http://www.clsi.org)

[standard@clsi.org](mailto:standard@clsi.org)

---

# One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test

Donna D. Castellone, MS, MASCP, MT(ASCP)SH  
Karen A. Moser, MD  
Elizabeth M. Van Cott, MD  
Richard A. Marlar, PhD  
Janet Cook, MT(ASCP)

Marilyn Johnston, ART  
Stephen Kitchen, FIBMS, PhD  
Samuel J. Machin, MB, ChB, FRCPath  
Diane Shafer, MT(ASCP)  
Laura Worfolk, PhD

---

## Abstract

Clinical and Laboratory Standards Institute document H47—*One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test* describes the principles and procedures necessary for the routine performance of the PT and APTT by conventional techniques using citrated plasma. Each of the two tests measures the time for a fibrin clot to develop in test plasma after activation. The chemical reactions are complex and, characteristically, results are affected by preexamination (preanalytical) and examination (analytical) variables. The PT and APTT are important screening tests used in laboratory evaluation of patients suspected to have disorders of blood coagulation, including the presence of circulating coagulation inhibitors. The PT measures the extrinsic or tissue factor pathway of the coagulation system and is used to monitor vitamin K antagonist therapy. The APTT measures the intrinsic coagulation pathway and is used in monitoring heparin therapy, as well as parenteral direct thrombin inhibitor anticoagulant therapy such as argatroban and lepirudin. The objective of this guideline is to improve test reproducibility through standardization of technique and ensure clinical relevance by setting test performance goals. The document also highlights the international effort for standardization of the PT through the use of the international normalized ratio.

Clinical and Laboratory Standards Institute (CLSI). *One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test*. 3rd ed. CLSI guideline H47 (ISBN 978-1-68440-172-7 [Print]; ISBN 978-1-68440-173-4 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2023.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org).

**If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:**

**P:** +1.610.688.0100 **F:** +1.610.688.0700 **E:** [customerservice@cls.org](mailto:customerservice@cls.org) **W:** [www.clsi.org](http://www.clsi.org)

Copyright ©2023 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, companion product, or other material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedure manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

## Suggested Citation

CLSI. *One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test*. 3rd ed. CLSI guideline H47. Clinical and Laboratory Standards Institute; 2023.

### Previous Editions:

June 1996, May 2008

H47-Ed3

ISBN 978-1-68440-172-7 (Print)

ISBN 978-1-68440-173-4 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 43, Number 5

.....

# Committee Membership

## Consensus Council

The Consensus Council sets priorities for CLSI standards development and votes on Final Draft documents to confirm that process requirements have been met. Consensus Council members are listed on the CLSI website: <https://clsi.org/standards-development/consensus-council/>

## Working Group on One-Stage PT/APTT

<p><b>Richard A. Marlar, PhD</b>  <b>Chairholder</b>  <b>Oklahoma City VA Medical Center</b>  <b>USA</b></p>	<p>Stephen Kitchen, FIBMS, PhD          Royal Hallamshire Hospital          United Kingdom</p>	<p>Diane Shafer, MT(ASCP)          Siemens Healthcare Diagnostics          USA</p>
<p>Janet Cook, MT(ASCP)          Coulter Diagnostics          USA</p>	<p>Samuel J. Machin, MB, ChB,          FRCPATH          The University College London          Hospitals          United Kingdom</p>	<p>Laura Worfolk, PhD          Quest Diagnostics Nichols Institute          USA</p>
<p>Marilyn Johnston, ART          Hemostasis Reference Laboratory          Canada</p>		

## Expert Panel on Hematology and Immunology and Ligand Assay

Expert panel volunteers support the development of CLSI documents by providing technical expertise in specialty areas. Expert panel members are listed by area of expertise on the CLSI website: <https://clsi.org/standards-development/expert-panels/>

## Staff

<p>Clinical and Laboratory Standards          Institute          USA</p>	<p>Laura Martin  <i>Editorial Manager</i></p>	<p>Kristy L. Leirer, MS  <i>Editor</i></p>
<p>Tabitha Kern, MS, MLS(ASCP)<sup>CM</sup>  <i>Program Manager</i></p>	<p>Catherine E.M. Jenkins, ELS  <i>Editor</i></p>	<p>Lisa M.W. Walker, MS, ELS  <i>Editor</i></p>

## Acknowledgment

---

CLSI, the Consensus Council, the Expert Panel on Hematology and Immunology and Ligand Assay, the Area Committee on Hematology, and the Working Group on One-Stage PT/APTT gratefully acknowledge the following volunteers for their important contributions to the limited revision of this guideline:

Donna D. Castellone, MS, MASCP,  
MT(ASCP)SH  
Columbia University  
USA

Karen A. Moser, MD  
ARUP Laboratories  
USA

Elizabeth M. Van Cott, MD  
Massachusetts General Hospital  
USA

### Acknowledgment in Memoriam of our Expert Panel on Hematology and Immunology and Ligand Assay Contributor and Colleague

The Expert panel on Hematology and Immunology and Ligand Assay would like to acknowledge the endless and valuable contributions made to this document by Dr. Elizabeth M. Van Cott. Dr. Van Cott was a true asset in the field of coagulation and laboratory medicine. Her passion for excellence and patient care was evident in all aspects of her work, and her contributions to this document will be a valuable tool for laboratories.

# Contents

.....

- Abstract ..... i
- Committee Membership..... iii
- Foreword..... ix
- Chapter 1: Introduction ..... 1
  - 1.1 Scope..... 2
  - 1.2 Standard Precautions..... 2
  - 1.3 Terminology..... 3
- Chapter 2: Equipment ..... 7
  - 2.1 Containers ..... 8
  - 2.2 Delivery Systems ..... 8
  - 2.3 Water Bath ..... 8
- Chapter 3: Specimen Collection, Transport, and Sample Storage ..... 9
- Chapter 4: Considerations in Performing the PT and APTT Tests ..... 11
  - 4.1 Manufacturer’s Instructions..... 12
  - 4.2 Acceptable Variability ..... 12
  - 4.3 Reagent Grade Water..... 12
  - 4.4 Calcium Ion Concentration ..... 12
  - 4.5 Conditions of the Test System ..... 12
  - 4.6 Controls Outside Stated Limits ..... 12
  - 4.7 Control Plasma Collection, Handling, and Storage..... 12
  - 4.8 Frequency of Control Testing ..... 13
  - 4.9 Reference Intervals ..... 13
  - 4.10 General Quality Control ..... 13
- Chapter 5: Performance of the PT Test ..... 15
  - 5.1 PT Principle ..... 16
  - 5.2 PT Reagents: Thromboplastins..... 16
  - 5.3 PT Performance Temperature..... 16
  - 5.4 PT Test Procedure..... 16
  - 5.5 PT End Point..... 16
  - 5.6 Validating and Calibrating PT Reagents ..... 17

.....

## Contents (Continued)

<b>Chapter 6: Local System Calibration</b>	<b>19</b>
6.1 General Guidelines for the Use of Certified Plasmas	20
6.2 Determination of Local ISI	20
6.3 Determination of Direct INR	21
6.4 Interpretation of Patient Test Results	22
6.5 PT Mixing Studies	23
6.6 Reporting PT Results	23
<b>Chapter 7: Performance of the APTT Test</b>	<b>25</b>
7.1 APTT Principle	26
7.2 APTT Reagent	26
7.3 APTT Performance Temperature	26
7.4 Contact Activation Time	26
7.5 APTT Test Procedure	26
7.6 APTT End Point	26
7.7 Unfractionated (UF) Heparin Sensitivity	27
7.8 Factor Sensitivity	27
7.9 Lupus Anticoagulants (LAs)	27
7.10 APTT Mixing Studies	28
7.11 Reporting of APTT Results	29
<b>Chapter 8: Monitoring Direct Thrombin Inhibitors and Other Anticoagulants</b>	<b>31</b>
8.1 Monitoring Parenteral Direct Thrombin Inhibitors	32
8.2 Measuring Direct Oral Anticoagulants	32
<b>Chapter 9: Effect of Hemophilia Therapies on APTT</b>	<b>33</b>
<b>Chapter 10: Sources of Error</b>	<b>35</b>
10.1 Specimen- or Sample-Related Problems	36
10.2 Possible Preexamination Errors	36
10.3 Possible Reagent-Related Problems	36
10.4 Possible Examination and Postexamination Errors	37

# Contents (Continued)

---

<b>Chapter 11: Supplemental Information</b> .....	<b>39</b>
<b>References</b> .....	<b>40</b>
<b>Appendix A.</b> Explanation of the International Sensitivity Index (ISI) .....	<b>46</b>
<b>Appendix B.</b> Description of the Geometric Mean (GM) .....	<b>47</b>
<b>Appendix C.</b> Establishing Heparin Therapeutic Range .....	<b>48</b>
<b>Appendix D.</b> Factor Sensitivity Determination .....	<b>50</b>
<b>The Quality Management System Approach</b> .....	<b>52</b>



## Foreword

---

Since its original description by Quick<sup>1</sup> in 1935, the prothrombin time (PT) has remained an important screening test in the laboratory evaluation of patients with suspected disorders of blood coagulation. It is the most common coagulation test performed in the clinical laboratory. Although the PT was originally described as a specific, one-stage assay of prothrombin or Factor II, it is sensitive to quantitative or qualitative abnormalities of any of the factors involved in the extrinsic and common pathways of the coagulation system (Factors II, V, VII, X, and fibrinogen), as well as inhibitors of these factors. It is an indicator of moderate to severe hepatic disease or chronic hepatic disease, as well as vitamin K deficiency or disseminated intravascular coagulation. The PT is also the most commonly used test for monitoring vitamin K antagonist (VKA) therapy. The PT may also be prolonged by direct oral anticoagulant (DOAC) therapy, depending on the nature and concentration of the DOAC used.

Thromboplastin, a phospholipid/tissue factor preparation and the principal reagent used in the PT assay, is commercially available in a variety of preparations of human or animal origin, or human or animal recombinant material. There are differences among commercial thromboplastin preparations in their responsiveness to reductions in coagulation factors that may affect their usefulness, particularly in the monitoring of VKA therapy.<sup>2-6</sup>

The activated partial thromboplastin time (APTT) is sensitive to quantitative and qualitative abnormalities in the intrinsic and common pathways of coagulation. It is the second most common coagulation procedure performed in routine laboratories. The APTT is particularly sensitive to defects of the intrinsic coagulation pathway (Factors VIII, IX, XI, XII, prekallikrein, and high-molecular-weight kininogen).<sup>7,8</sup> It is commonly used for monitoring unfractionated heparin anticoagulant therapy. It detects other types of pathological inhibitors of blood coagulation, the most common of which is the lupus anticoagulant (LA), and it is used to monitor factor replacement therapy. APTT reagents are a mixture of procoagulant phospholipids and a contact activator. The phospholipids may be of human, animal, or vegetable origin, and there are a variety of activating substances (eg, celite, kaolin, micronized silica, ellagic acid).

Ideally, the APTT is prolonged when levels of coagulation factor activity are sufficiently decreased to cause clinical bleeding (<30%).<sup>9</sup> However, a number of studies have shown considerable differences in the responsiveness of the various APTT reagents to mild and moderate factor deficiencies, particularly deficiencies of Factor VIII and/or Factor IX.<sup>7,8,10,11</sup> A similarly variable sensitivity of the APTT to circulating LAs has been reported.<sup>12</sup> Likewise, marked APTT variability in responsiveness to heparin and DOAC therapy has been observed among commercially available APTT reagents.<sup>8,13,14</sup>

This document is written for laboratory and/or clinical personnel responsible for the performance, quality control, and reporting of the PT and APTT tests, as well as for manufacturers of coagulation instruments and reagents who are responsible for maintaining appropriate performance standards. This document should be used in conjunction with CLSI documents H54<sup>15</sup> and H57.<sup>16</sup>

## Overview of Changes

This guideline was revised in 2023 under the Limited Revision Process and replaces the second edition of the guideline, which was published in 2008. Several changes were made in this edition, including:

- Adding new anticoagulants
- Adding effect of hemophilia therapies on APTT
- Adding new references

**NOTE:** The content of this guideline is supported by the CLSI consensus process and does not necessarily reflect the views of any single individual or organization.

### KEY WORDS

activated partial  
thromboplastin time (APTT)

citrate

coagulation

coagulation factor(s)

control (plasma)

fibrinogen

international normalized ratio  
(INR)

international sensitivity index  
(ISI)

phospholipids

prothrombin time (PT)

thrombin time

thromboplastin

tissue factor



# Chapter ①

## Introduction

# One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test

---

## 1 Introduction

The results of the PT and APTT tests can be affected by a number of preexamination variables, such as method of blood collection; surface characteristics of collection containers; type and concentration of anticoagulant; specimen and sample storage conditions; and examination variables, such as sample incubation time and temperature, contact activation time, type of reagents, and the method of end-point detection. In this document, standard methods for collection, transport, and processing of blood specimens are referenced in CLSI document H21,<sup>17</sup> and test performance specifications are described. This is intended to minimize the effects of such variables, improve precision and accuracy, and, thus, the clinical usefulness of the PT and APTT.

### 1.1 Scope

H47 provides general guidelines for performing PT and APTT tests with the conventional method of using citrated, platelet-poor plasma. It does not cover alternative methods of using citrated whole blood or capillary blood obtained by the fingerstick method. Nonclotting-based, end-point-detection tests, such as chromogenic substrate assays, are also not discussed.

### 1.2 Standard Precautions

Because it is often impossible to know what isolates or specimens might be infectious, all patient and laboratory specimens are treated as infectious and handled according to “standard precautions.” Standard precautions are guidelines that combine the major features of “universal precautions and body substance isolation” practices. Standard precautions cover the transmission of all infectious agents and thus are more comprehensive than universal precautions, which are intended to apply only to transmission of blood-borne pathogens. Standard and universal precaution guidelines are available from the US Centers for Disease Control and Prevention.<sup>18</sup> For specific precautions for preventing the laboratory transmission of all infectious agents from laboratory instruments and materials and for recommendations for the management of exposure to all infectious disease, refer to CLSI document M29.<sup>19</sup>