



CLINICAL AND  
LABORATORY  
STANDARDS  
INSTITUTE®

2nd Edition

# M24S

## Performance Standards for Susceptibility Testing of *Mycobacteria*, *Nocardia* spp., and Other Aerobic Actinomycetes

This document includes updated breakpoint and quality control tables for the Clinical and Laboratory Standards Institute susceptibility testing standard M24.

A CLSI supplement for global application.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advances in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeal Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeal, documented in the CLSI *Standards Development Policies and Processes*, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute

P: +1.610.688.0100

F: +1.610.688.0700

[www.clsi.org](http://www.clsi.org)

[standard@clsi.org](mailto:standard@clsi.org)

# Performance Standards for Susceptibility Testing of *Mycobacteria, Nocardia spp.*, and Other Aerobic Actinomycetes

Nicole M. Parrish, PhD, MHS, D(ABMM)  
Nancy L. Wengenack, PhD, D(ABMM)  
Adam Barker, PhD  
Barbara A. Brown-Elliott, MS, MT(ASCP)SM  
Daniela Maria Cirillo, MD, PhD  
Susan Harrington, PhD, D(ABMM), MLS(ASCP)<sup>CM</sup>  
Reeti Khare, PhD, D(ABMM)  
Scott B. Killian, BS  
Richard Pfeltz, PhD  
Elvira Richter, PhD  
Marie-Claire Rowlinson, PhD, D(ABMM)  
Adrian M. Zelazny, PhD, D(ABMM)

## Abstract

Clinical and Laboratory Standards Institute document M24S—*Performance Standards for Susceptibility Testing of Mycobacteria, Nocardia spp., and Other Aerobic Actinomycetes* includes the minimal inhibitory concentrations and QC ranges developed following the standards described in CLSI document M24.<sup>1</sup> The data in the tables are valid only when the methodology in CLSI document M24<sup>1</sup> is followed.

Clinical and Laboratory Standards Institute (CLSI). *Performance Standards for Susceptibility Testing of Mycobacteria, Nocardia spp., and Other Aerobic Actinomycetes*. 2nd edition. CLSI supplement M24S (ISBN 978-1-68440-168-0 [Print]; ISBN 978-1-68440-169-7 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2023.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org).

**If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:**

**P:** +1.610.688.0100 **F:** +1.610.688.0700 **E:** [customerservice@clsi.org](mailto:customerservice@clsi.org) **W:** [www.clsi.org](http://www.clsi.org)

Copyright ©2023 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, derivative product, or other material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedures manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

## Suggested Citation

CLSI. *Performance Standards for Susceptibility Testing of Mycobacteria, Nocardia spp., and Other Aerobic Actinomycetes*. 2nd edition. CLSI supplement M24S. Clinical and Laboratory Standards Institute; 2023.

## Previous Editions:

M62: November 2018

M24-Ed3-S-Ed2

ISBN 978-1-68440-168-0 (Print)

ISBN 978-1-68440-169-7 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 43, Number 2

## Committee Membership

### Working Group on Antimycobacterial Susceptibility Testing

Nicole M. Parrish, PhD, MHS,  
D(ABMM)  
Co-Chairholder  
Johns Hopkins Hospital-Pathology  
USA

Nancy L. Wengenack,  
PhD, D(ABMM)  
Co-Chairholder  
Mayo Clinic  
USA

Adam Barker, PhD  
ARUP Laboratories  
USA

Barbara A. Brown-Elliott,  
MS, MT(ASCP)SM  
University of Texas Health  
Science Center at Tyler  
USA

Daniela Maria Cirillo, MD, PhD  
World Health Organization  
Collaborating Centre and TB  
Supranational Reference Laboratory  
Italy

Susan Harrington, PhD, D(ABMM),  
MLS(ASCP)<sup>CM</sup>  
Cleveland Clinic  
USA

Reeti Khare, PhD, D(ABMM)  
National Jewish Health  
USA

Scott B. Killian, BS  
Thermo Fisher Scientific  
USA

Richard Pfeltz, PhD  
Becton Dickinson and Company  
USA

Elvira Richter, PhD  
Laboratory Limbach  
Germany

Marie-Claire Rowlinson,  
PhD, D(ABMM)  
New York State Department  
of Health  
USA

Adrian M. Zelazny, PhD, D(ABMM)  
National Institutes of Health  
Department of Laboratory Medicine  
USA

### Staff

Clinical and Laboratory Standards  
Institute  
USA

Christine M. Lam, MLS(ASCP)  
*Program Manager*

Laura Martin  
*Editorial Manager*

Catherine E.M. Jenkins, ELS  
*Editor*

Kristy L. Leirer, MS  
*Editor*

Lisa M.W. Walker, MS, ELS  
*Editor*

### Acknowledgment

CLSI and the Working Group on Antimycobacterial Susceptibility Testing gratefully acknowledge the following volunteers for their important contributions to the development of this document:

Lynette Y. Berkeley, PhD, MT(ASCP)  
FDA Center for Drug Evaluation and  
Research  
USA

Patricia S. Conville, MS, MT  
FDA Center for Devices and  
Radiological Health  
USA

Chris Lewis  
Thermo Fisher Scientific  
USA

Chris Pillar, PhD  
Microbiologics  
USA

Max Salfinger, MD, FAAM,  
FIDSA, ATSF  
University of South Florida  
College of Public Health & Morsani  
College of Medicine  
USA



## Contents

Abstract .....	i
Committee Membership.....	iii
Overview of Changes .....	vii
Abbreviations and Acronyms .....	xi
Reference .....	xii
Table 1. Broth Microdilution Breakpoints and Interpretive Categories for MTBC Tested in Middlebrook 7H9 Broth Supplemented With OADC Using MIC Panels.....	1
Table 2. Susceptibility Testing of MTBC to Second-Line Drugs Using the Fluorescence-based Commercial Shorter-Incubation Liquid Media System.....	3
Table 3. Antimicrobial Agents and Breakpoints for Testing MAC.....	6
Table 4. Antimicrobial Agents and Breakpoints for Testing <i>Mycobacterium kansasii</i> .....	8
Table 5. Antimicrobial Agents and Breakpoints for Testing Slowly Growing Nontuberculous <i>Mycobacteria</i> Other Than MAC and <i>Mycobacterium kansasii</i> .....	9
Table 6. Antimicrobial Agents and Breakpoints for Testing Rapidly Growing <i>Mycobacteria</i> .....	10
Table 7. Antimicrobial Agents and Breakpoints for Testing <i>Nocardia</i> spp. and Other Aerobic Actinomycetes.....	13
Table 8. Expected Antimicrobial Susceptibility Patterns of the Most Commonly Isolated <i>Nocardia</i> spp. ....	15
Table 9. Expected Results for Routine and Supplemental QC Strains When Testing MTBC by a Culture-based Method.....	17
Table 10. MIC QC Ranges When Testing MTBC Using <i>Mycobacterium tuberculosis</i> ATCC® 27294 (H37Rv) Tested in Middlebrook 7H9 Medium Supplemented With OADC .....	18
Table 11. QC Ranges for Susceptibility Testing of <i>Mycobacterium tuberculosis</i> ATCC® 27294 (H37Rv) to Second-Line Drugs Using the Fluorescence-based Commercial Shorter-Incubation Liquid Media System .....	20
Table 12. Broth Microdilution QC Ranges When Testing Slowly Growing Nontuberculous <i>Mycobacteria</i> in Cation-Adjusted Mueller-Hinton Broth With OADC .....	21
Table 13. Broth Microdilution QC Ranges When Testing Rapidly Growing <i>Mycobacteria</i> in Cation-Adjusted Mueller-Hinton Broth Without OADC .....	23

## Contents (Continued)

Table 14. Broth Microdilution QC Ranges When Testing <i>Nocardia</i> spp. and Other Aerobic Actinomycetes in Cation-Adjusted Mueller-Hinton Broth Without OADC .....	25
The Quality Management System Approach .....	27

## Overview of Changes

This document replaces the previous edition of the approved document, M62-Ed1, published in 2018. **This supplement contains the most current information and serves as a means to notify laboratories of important changes that occur between published editions of CLSI document M24.**<sup>1</sup> Several changes were made in this edition, including:

Section/Table	Action	Change to:	Reason/Specific Change
Key Words	Added	Key words	<i>Mycobacterium avium</i> complex
General	Revised	Footnotes	Footnotes were reorganized as general and specific comments.
Table 1. Broth Microdilution Breakpoints and Interpretive Categories for MTBC Tested in Middlebrook 7H9 Broth Supplemented With OADC Using MIC Panels	Added	Comments and associated references	<ul style="list-style-type: none"> <li>Regarding important changes between this document and CLSI document M24<sup>1</sup></li> <li>Regarding establishment of breakpoints using lyophilized MIC panels, their incubation requirements, and importance of following manufacturer's media formulations</li> <li>Regarding testing for rifampin susceptibility using automated broth-based systems</li> </ul>
	Revised	Comment (5)	NOTE regarding repeat testing clarified
		MIC breakpoints for rifampin	Updated WHO standard concentration
	Deleted	Title	The term "commercially available"
Table 2. Susceptibility Testing of MTBC to Second-Line Drugs Using the Fluorescence-based Commercial Shorter-Incubation Liquid Media System	Revised	Comment	Clarified amikacin and kanamycin guidance
Table 3. Antimicrobial Agents and Breakpoints for Testing MAC	Added	Alternate routine QC strains and comment	For clarithromycin, linezolid, and moxifloxacin

## Overview of Changes (Continued)

Section/Table	Action	Change to:	Reason/Specific Change
Table 4. Antimicrobial Agents and Breakpoints for Testing <i>Mycobacterium kansasii</i>	Added	Alternate routine QC strains and comment	For clarithromycin, linezolid, and moxifloxacin
Table 5. Antimicrobial Agents and Breakpoints for Testing Slowly Growing Nontuberculous Mycobacteria Other Than MAC and <i>M. kansasii</i>	Added	Alternate routine QC strains and comment	For clarithromycin, linezolid, and moxifloxacin
Table 6. Antimicrobial Agents and Breakpoints for Testing Rapidly Growing Mycobacteria	Revised	Comment	<ul style="list-style-type: none"> <li>• Included newly described species: <ul style="list-style-type: none"> <li>– <i>Mycobacterium conceptionense</i></li> <li>– <i>Mycobacterium farcinogenes</i></li> <li>– <i>Mycobacterium syngnathidarum</i></li> </ul> </li> </ul>
Table 8. Expected Antimicrobial Susceptibility Patterns of the Most Commonly Isolated <i>Nocardia</i> spp.	Added	Comments	<ul style="list-style-type: none"> <li>• General comment regarding <i>Nocardia</i> susceptibility patterns</li> <li>• Interpretive comment for <i>Nocardia brasiliensis</i></li> </ul>
		References	Regarding susceptibility testing and profiling of <i>Nocardia</i> spp.
	Revised	Expected patterns	<ul style="list-style-type: none"> <li>• Various agents with: <ul style="list-style-type: none"> <li>– <i>N. brasiliensis</i></li> <li>– <i>Nocardia cyriacigeorgica</i></li> <li>– <i>Nocardia farcinica</i></li> <li>– <i>Nocardia nova</i> complex</li> <li>– <i>Nocardia otitidiscaviarum</i></li> <li>– <i>Nocardia pseudobrasiliensis</i></li> </ul> </li> </ul>
Table 10. MIC QC Ranges When Testing MTBC Using <i>Mycobacterium tuberculosis</i> ATCC <sup>®a</sup> 27294 (H37Rv) Tested in Middlebrook 7H9 Medium Supplemented With OADC	Added	Comment	Regarding the establishment of breakpoints using lyophilized or frozen MIC panels, incubation requirements, and importance of following manufacturer's media formulations
	Deleted	General	The term "commercially available"

## Overview of Changes (Continued)

Section/Table	Action	Change to:	Reason/Specific Change
Table 12. Broth Microdilution QC Ranges When Testing Slowly Growing Nontuberculous Mycobacteria in Cation-Adjusted Mueller-Hinton Broth With OADC	Added	Title	Text specifying broth used
		MIC QC ranges, comments, and associated reference	Fobrepodacin tested with routine and supplemental QC strains
		Comment	Regarding alternate QC strains for clarithromycin, linezolid, and moxifloxacin
		Comment	Regarding morphological variations observed with <i>M. avium</i> ATCC® 700898
	QC strain	<i>M. avium</i> ATCC® 700898 as an alternate QC organism for MAC testing for amikacin, clarithromycin, linezolid, and moxifloxacin	
Revised	MIC QC ranges	<ul style="list-style-type: none"> <li>• <i>Staphylococcus aureus</i> ATCC® 29213 <ul style="list-style-type: none"> <li>– Moxifloxacin</li> <li>– Rifampin</li> </ul> </li> </ul>	
Table 13. Broth Microdilution QC Ranges When Testing Rapidly Growing Mycobacteria in Cation-Adjusted Mueller-Hinton Broth Without OADC	Added	Title	Text specifying broth used
		MIC QC ranges, comments, and associated reference	Fobrepodacin tested with routine and supplemental QC strains
	MIC QC ranges	Omadacycline tested with routine and supplemental QC strains	
Revised	MIC QC ranges	<ul style="list-style-type: none"> <li>• <i>Pseudomonas aeruginosa</i> ATCC® 27853 <ul style="list-style-type: none"> <li>– Ciprofloxacin</li> <li>– Meropenem</li> </ul> </li> <li>• <i>S. aureus</i> ATCC® 29213 <ul style="list-style-type: none"> <li>– Imipenem</li> </ul> </li> </ul>	
Table 14. Broth Microdilution QC Ranges When Testing <i>Nocardia</i> spp. and Other Aerobic Actinomycetes in Cation-Adjusted Mueller-Hinton Broth Without OADC	Added	Title	Text specifying broth used
	Revised	MIC QC range	<i>S. aureus</i> ATCC® 29213 and imipenem

Abbreviations. ATCC®, American Type Culture Collection; MAC, *Mycobacterium avium* complex; MIC, minimal inhibitory concentration; MTBC, *Mycobacterium tuberculosis* complex; OADC, oleic acid-albumin-dextrose-catalase; QC, quality control; WHO, World Health Organization.

## Overview of Changes (Continued)

### Footnote

- a. ATCC® is a trademark of the American Type Culture Collection.

**NOTE:** The content of this document is supported by the CLSI consensus process and does not necessarily reflect the views of any single individual or organization.

### Key Words

aerobic actinomycetes, antimicrobial agents, antimicrobial susceptibility testing, antituberculous drugs, ***Mycobacterium avium complex***, *Mycobacterium tuberculosis* complex, *Nocardia* spp., nontuberculous mycobacteria

.....

## Abbreviations and Acronyms

<b>ATCC<sup>®a</sup></b>	American Type Culture Collection
<b>CO<sub>2</sub></b>	carbon dioxide
<b>HPLC</b>	high-performance liquid chromatography
<b>I</b>	intermediate
<b>IV</b>	intravenous
<b>MAC</b>	<i>Mycobacterium avium</i> complex
<b>MIC</b>	minimal inhibitory concentration
<b>MTBC</b>	<i>Mycobacterium tuberculosis</i> complex
<b>OADC</b>	oleic acid-albumin-dextrose-catalase
<b>QC</b>	quality control
<b>QMS</b>	quality management system
<b>QSE</b>	quality system essential
<b>R</b>	resistant
<b>S</b>	susceptible
<b>V</b>	variable
<b>WHO</b>	World Health Organization

---

<sup>a</sup> ATCC<sup>®</sup> is a registered trademark of the American Type Culture Collection.

.....

## Reference

- <sup>1</sup> CLSI. *Susceptibility Testing of Mycobacteria, Nocardia spp., and Other Aerobic Actinomycetes*. 3rd ed. CLSI standard M24. Clinical and Laboratory Standards Institute; 2018.

**Table 1. Broth Microdilution Breakpoints and Interpretive Categories for MTBC Tested in Middlebrook 7H9 Broth Supplemented With OADC Using MIC Panels<sup>1-3</sup>**

QC Recommendation (see Table 10 for acceptable QC ranges)

Routine QC strain:

- *Mycobacterium tuberculosis* ATCC® 27294 (H37Rv)

#### General Comments

- (1) ATCC® is a registered trademark of the American Type Culture Collection.
- (2) This table contains updated information and notifies laboratories of important changes among published versions of CLSI document M24.<sup>4</sup>
- (3) The breakpoints in Table 1 were established using lyophilized MIC panels incubated in 5% to 10% CO<sub>2</sub>. Incubating tests in ambient air may result in erroneous MICs because MTBC grow more slowly in ambient air. If preparing in-house MIC panels, care should be taken to match the manufacturer's formulation because the use of polysorbate or glycerol in the broth may affect the MIC values obtained.
- (4) This antimicrobial susceptibility testing system is not regulatory organization cleared.

**NOTE: Information in black boldface type is new or modified since the previous edition.**

Antimicrobial Agent	MIC Breakpoints, µg/mL			Comments
	S	Inconclusive <sup>3</sup>	R	
Ethambutol	≤ 2	4	≥ 8	(5) Inconclusive MIC for ethambutol. An MIC of 4 µg/mL obtained by broth microdilution using <b>lyophilized panels</b> does not correlate with either a susceptible or resistant result in commercial automated, short-incubation broth systems, and there are no clinical data correlating such a result with ethambutol treatment response. <b>NOTE: Repeat testing using an alternative broth method (eg, critical concentration) or genotypic method</b> may determine whether the isolate in question is susceptible or resistant.